

Sample S001 | LKM, --

Methodics	Negative	Total
ALEGRIA, Orgentec	2	2
BlueDiver Immunodot, D-tek	1	1
EliA Phadia, Thermo Fisher	3	3
ELISA, Diesse Chorus	1	1
ELISA, Euroimmun	3	3
ELISA, Inova	1	1
Euroline, Euroimmun	21	21
IFA, BioSystems	1	1
IFA, Euroimmun	15	15
IFA, Inova	3	3
IFA Kallestadt, Bio-Rad	1	1
IFA Zenit, Menarini	1	1
Microblot-Array Liver-Profile,TestLine	1	1
Total	54	54

Sample S001 | SMA, --

Methodics	Negative	Weak positive	Positive	Total
ELISA, Euroimmun	1	-	-	1
ELISA, Inova	1	-	-	1
Euroline, Euroimmun	4	-	-	4
IFA, AESKU DIAGNOSTICS AESKUSLIDES rLKS	1	-	-	1
IFA, BioSystems	-	-	1	1
IFA, Euroimmun	26	-	-	26
IFA, GA Generic Assays	1	-	-	1
IFA, Inova	5	1	-	6
IFA Kallestadt, Bio-Rad	1	-	-	1
IFA Zenit, Menarini	1	-	-	1
Total	41	1	1	43

Sample S001 | AMA, --

Methodics	Negative	Total
ALEGRIA, Orgentec	6	6
BlueDiver Immunodot, D-tek	1	1
ELiA Phadia, Thermo Fisher	4	4
ELISA, Dienes Chorus	1	1
ELISA, Euroimmun	2	2
ELISA, Inova	1	1
Euroline, Euroimmun	15	15
IFA, AESKU DIAGNOSTICS AeskusLIDES rLKS	1	1
IFA, BioSystems	1	1
IFA, Euroimmun	21	21
IFA, Inova	4	4
IFA Kallestadt, Bio-Rad	1	1
IFA Zenit, Menarini	1	1
Microblot-Array Liver-Profile,TestLine	1	1
QuantaLite M2 EP (MIT3) ELISA	1	1
Snibe Maglumi AMA-M2 IgG (CLIA)	1	1
Total	62	62

Sample S001 | PCA, --

Methodics	Negative	Positive	Total
ALEGRIA, Orgentec	4	-	4
BlueDiver Immunodot, D-tek	1	-	1
EliA Phadia, Thermo Fisher	9	-	9
ELISA, Euroimmun	2	-	2
ELISA, Inova	1	-	1
ELISA, Parietalcell-Ab, IBL	2	-	2
IFA, AESKU DIAGNOSTICS AESKUSLIDES rLKS	2	-	2
IFA, BioSystems	1	-	1
IFA, Euroimmun	17	1	18
IFA, Inova	2	-	2
IFA Kallestadt, Bio-Rad	1	-	1
IFA Zenit, Menarini	1	-	1
Microblot-Array Liver-Profile,TestLine	-	1	1
Total	43	2	45

Sample S002 | LKM, --

Methodics	Negative	Total
ALEGRIA, Orgentec	2	2
BlueDiver Immunodot, D-tek	1	1
EliA Phadia, Thermo Fisher	3	3
ELISA, DIESSE Chorus	1	1
ELISA, Euroimmun	3	3
ELISA, Inova	1	1
Euroline, Euroimmun	21	21
IFA, BioSystems	1	1
IFA, Euroimmun	15	15
IFA, Inova	3	3
IFA Kallestadt, Bio-Rad	1	1
IFA Zenit, Menarini	1	1
Microblot-Array Liver-Profile,TestLine	1	1
Total	54	54

Sample S002 | SMA, --

Methodics	Negative	Positive	Total
ELISA, Euroimmun	1	-	1
ELISA, Inova	-	1	1
Euroline, Euroimmun	2	2	4
IFA, AESKU DIAGNOSTICS AESKUSLIDES rLKS	-	1	1
IFA, BioSystems	-	1	1
IFA, Euroimmun	13	13	26
IFA, GA Generic Assays	-	1	1
IFA, Inova	3	3	6
IFA Kallestadt, Bio-Rad	-	1	1
IFA Zenit, Menarini	-	1	1
Total	19	24	43

Sample S002 | AMA, --

Methodics	Negative	Total
ALEGRIA, Orgentec	6	6
BlueDiver Immunodot, D-tek	1	1
EliA Phadia, Thermo Fisher	4	4
ELISA, DIESSE Chorus	1	1
ELISA, Euroimmun	2	2
ELISA, Inova	1	1
Euroline, Euroimmun	15	15
IFA, AESKU DIAGNOSTICS AESKUSLIDES rLKS	1	1
IFA, BioSystems	1	1
IFA, Euroimmun	21	21
IFA, Inova	4	4
IFA Kallestadt, Bio-Rad	1	1
IFA Zenit, Menarini	1	1
Microblot-Array Liver-Profile,TestLine	1	1
QuantaLite M2 EP (MIT3) ELISA	1	1
Snibe Maglumi AMA-M2 IgG (CLIA)	1	1
Total	62	62

Sample S002 | PCA, --

Methodics	Negative	Positive	Total
ALEGRIA, Orgentec	4	-	4
BlueDiver Immunodot, D-tek	1	-	1
EliA Phadia, Thermo Fisher	9	-	9
ELISA, Euroimmun	2	-	2
ELISA, Inova	1	-	1
ELISA, Parietalcell-Ab, IBL	2	-	2
IFA, AESKU DIAGNOSTICS AESKUSLIDES rLKS	2	-	2
IFA, BioSystems	1	-	1
IFA, Euroimmun	16	2	18
IFA, Inova	2	-	2
IFA Kallestadt, Bio-Rad	1	-	1
IFA Zenit, Menarini	1	-	1
Microblot-Array Liver-Profile,TestLine	1	-	1
Total	43	2	45

Report info

Participants

57 participants from 16 countries.

Report info

The results are divided into groups according to the method stated by the laboratory and presented in laboratory-specific tables. Accepted results are marked with green colour and laboratory’s own result with a black radio button . If you have not reported any results you will get a note: “You have not responded in time, only global report is available.” Always compare your results to those obtained with the same test method.

For information on report interpretation and performance evaluation, please see the "EQAS Interpretation guidelines" in LabScala User instructions. In case you have any questions regarding the reports, please contact the EOA Coordinator.